**BLANCONE<sup>®</sup>** **ARCUS<sup>+</sup>**

**FIXED**

**LED Bleaching Lamp**

**Lampada sbiancante a LED**

**LED Bleaching Lampe**

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## LED Bleaching lamp FIXED

GB

<b>I.</b>	Unit description and functions .....	2
<b>II.</b>	Technical specifications.....	2
<b>III.</b>	Packaging / components .....	3
<b>IV.</b>	Safety precautions .....	3
<b>V.</b>	Installation.....	6
<b>VI.</b>	Preparation for operation .....	8
<b>VII.</b>	Daily care and maintenance .....	9
<b>VIII.</b>	Preparation for work with patients.....	12
<b>IX.</b>	Positions of bleaching head .....	12
	against patient's teeth	
<b>X.</b>	Determination of distance .....	13
	to patient's teeth	
<b>XI.</b>	Problems and solutions .....	13
<b>XII.</b>	Warranty .....	14
<b>XIII.</b>	Symbols .....	85
<b>XIV.</b>	Service data .....	85
<b>XV.</b>	Declaration of conformity .....	86
<b>XVI.</b>	Led bleaching unit data .....	87

Caution!

Before installing and operating with the unit, carefully read this manual.

## I. UNIT DESCRIPTION AND FUNCTIONS

BLANCONE ARCUS bleaching device is designed to bleach the solid tooth tissue by photoactivation of a whitening gel applied in dentistry. The appliance should only be used by a dental practitioner and in a dental practice. The device consists of a Bleaching Head, Stand, Fixing bracket and a Power Adapter. BLANCONE ARCUS BUILT-IN is manufactured in accordance with the requirements of Directive 93/42/EEC.

## II. TECHNICAL SPECIFICATIONS

1. Operating voltage to the power adapter - 100-240V / 50-60Hz, to the whitening head - 24VDC.
2. Current consumption - 0.4-0.2 A
3. Dimensions:
  - Bleaching head
    - length - 220 mm
    - width - 155 mm
    - height - 95 mm
  - Stand
    - length (assembled) - 52 cm / 80 cm with extended telescopic arm
    - height 26 cm
    - maximum diameter of the standard fixing bracket - 50 mm.
    - Optional fixing brackets with max. diameters 60mm, 70mm, 80mm and 90mm
4. Weight:
  - Bleaching head - 640 g
  - Stand - 2000 g
5. Irradiance - over 250 mW / sq.cm. measured in the output window.
6. Operation time - 10 sec. - 30 min /  $\pm$  5% /
7. Emitted light - blue 430 - 490 nm.
8. Light source - 10 LED modules x 5W = 50W
9. Fan air cooling.
10. Ability to start and stop the Bleaching Head at any time.
11. Maximum continuous operation time 99 min – after every 99 minutes of continuous operation let the unit cool for at least 10 minutes.

Upon request the manufacturer of this device will provide all other necessary documentation and/or information for the user's technical staff to repair those parts of the apparatus which the manufacturer has indicated as suitable for repair outside the manufacturer's aftermarket service facility.

### III. PACKAGING / COMPONENTS

1. Power cord 24 V DC to the Bleaching Head
2. Straight arm
3. Telescopic arm
4. Central position locking mechanism
5. Angular Arm
6. Fixing bracket
7. Bleaching head positioning mechanism
8. Bleaching head
9. Dental unit vertical arm
10. Control panel
- 10A. "+" Timer or "Pause" button
- 10B. "-" Timer or "Pause" button
- 10C. Start / Stop button
- 10D. Display
11. Cable retaining clip
12. Cables and connectors
13. Safety glasses - 2 pcs.
14. Power Adapter 100-240V AC / 24V DC.

### IV. SAFETY PRECAUTIONS

1. The device must be used in accordance with the User Manual. It is not allowed to use it for any other purpose or by unauthorized personnel.
2. The device is equipped with European type plug for voltage 100-240V/50-60Hz. Usage of the device with other plugs and/or other voltages can damage it or may lead to certain malfunctions.

Electrical safety is provided by Class I protection against electric shock according to EN 60601-1.

3. BLANCONE ARCUS should be operated indoors under the following conditions:
  - temperature from + 10 ° to + 40 ° C;
  - relative humidity of the air - 30 - 75%; The device should not be used in rooms with high humidity levels, condensation or wet floor, walls and/or ceiling.
  - atmospheric pressure 700 - 1060 hPa;
  - absence of chemically active and flammable substances.
4. Strong electromagnetic fields in the building may cause interference and malfunction. If their source cannot be determined, change the location and turn it on in another electrical plug.
5. BLANCONE ARCUS is a source of extremely intense light in the blue light spectrum to which the human eye has a high sensitivity. Serious measures need to be taken in order to protect patients, medical staff and accidentally close people, animals and plants. Always use protective glasses for the operator and for the patient's eyes, wear mask and facial sunscreen. Eye irritation carries a high risk. Never direct the light to the eyes! Irradiation must be limited to the place of work.

#### Caution!

The light emitted is at a wavelength of 430-490 nanometers.

Use safety glasses that comply with the following requirements:

- tightly cover the eyes and the temple even if the operator bears dioptric glasses.
  - made from a voluminous colored impact-resistant plastic.
  - admit light with a wavelength exceeding 600nm.
  - reduce the intensity of the blue spectrum over 100 times.
  - have a stable mechanical structure and do not have any scratches or damage on their surface.
6. The device should not be used with pregnant and breast-feeding mothers, patients with: oncological problems, severe parodontal pathology, recessions, dental paresthesia and under 17 years of age.  
The device may only be used after medical consultation on or by persons with implanted cardiac pacemaker; persons suffering from photobiological reactions; taking photosensitive medications; individuals with cataract surgery, people with retinal illnesses, allergy sufferers; persons with multiple sensitive skin or dermatitis, etc.  
If the bleaching protocol is not followed, pain, hypersensitivity, enamel defects and even burning of the non-calcified, soft tissue may occur.
  7. It is desirable to check the sensitivity of the patient to the intensity of the light. The device has two levels of intensity (High and Low) and the optimal conditions are selected by adjusting the distance of the device to the patient's teeth. Testing should begin at a high level "Hi" (High) and a distance of 5 cm. If the thermal sensation is unpleasant and irritating, increase the distance to achieve an acceptable thermal comfort. If the thermal sensation is still unpleasant at 10-15 cm switch to low intensity - "Lo" (Low) and shorten the

distance to the teeth. If the patient has no sensitivity or sensitivity cannot be measured for some reason, longer distance should be used.

8. Do not be placed any solvents, flammable liquids and heat sources near the device and/ or the power adapter, because exposing it to any of these may cause damage, fire or an explosion.
9. Cleaning agents should not penetrate the device, as this may lead to malfunctions, damage or fire.
10. Protect the outer insulation of all cables from damage that may be caused by sharp objects, strong pulling, rodents and chemical substances. If any such damages are noticed switch off the device and immediately take it to the service facility.
11. In the event of a thunderstorm, the device must be stopped and disconnected from the power network.
12. Do not force the horizontal arm to avoid mechanical shocks to the counterweight or the bleaching head. In the event of a mechanical shock, if the bleaching head is damaged, the device should not be used. It must be taken immediately to the service facility.
13. The movement of the bleaching head and horizontal arm without adjustments to the Central and/ or Front position locking mechanism may damage the device.
14. Any repairs may only be performed by service technicians authorized by the manufacturer.
15. When replacing damaged components, only original BLANCONE ARCUS spare parts should be used. The warranty of the device does not cover damages caused by the use of non-original spare parts. The device or any of its components should not be disassembled while plugged into the power network!
16. Before each procedure, the bleaching head window should be disinfected with disinfecting solution.
17. Take precautions against accidental injury to the human body from mechanical parts (moving and stationary), pinching, inertial acceleration of the counterweight etc. Safe handling of the device is vital to avoid injury.
18. Before the procedure each patient must sign an informed consent, detailing the possible side effects of the process. Note: It is also advisable to visit a suitable Teeth whitening in dental practice course.
19.  Fragile! Use caution when transporting, using and storing the device! When transported in an assembled condition, two persons are required to keep the moving parts. With any unevenness on the flooring, the device should only be moved by two people.
20.  According to Directive 2012/19/EEC, this symbol indicates that the product should not be disposed as a general waste at the end of its lifespan. The product must be taken to a specialized center for the separate collection of electrical and electronic equipment. Proper disposal of equipment that is no longer used prevents negative consequences for the environment and human health!

## V. INSTALLATION

BLANCONE ARCUS BUILT-IN is designed to be mounted with a fixing bracket to the vertical arm of the dental unit.

1. Unpack the transport box, remove the components and initially assemble the fixing bracket to the dental unit. Make sure the diameter of the bracket matches the diameter of the dental unit vertical arm. Soft PVC tape strips are provided in the case that there is a need to compensate for small difference in the two diameters. Tighten the two screws with a 4 mm Allen key (included in the package).

2. Tighten the stop/safety screw with a 3 mm Allen key (included in the package).

3. Insert the cable through the small opening of the fixing bracket and place the angular arm in it.


4. Attach the straight arm and the cable retaining clip. Tighten the stop/safety screw with a 3 mm Allen key (included in the package). Make sure that the straight arm can be moved up and down for precise positioning.

5. Altering the force required to move the straight arm up and down is done by adjusting the screw (located on the back of the central position locking mechanism) with a 5 mm Allen key (included in the package).

6. Mount the bleaching head by turning clockwise the positioning knob located on the telescopic arm and fix the cable with the retaining clip.


7. Connect the cable from the top of the angular arm to the one coming from the bleaching head using the connectors (see figure). Note: the device will not work if the cable connectors are incorrectly coupled.

8. Connect the cable on the bottom of the angular arm with the power adapter using the connectors (see figure). Note: the device will not work if the cable connectors are incorrectly coupled.


9. Plug the power cord into the power network.

## VI. PREPARATION FOR OPERATION

1. The BLANCONE ARCUS bleaching device should be placed in a suitable position around the patient and the dental unit. See the positioning examples (p.IX). The bleaching head must be in the same vertical plane and parallel with the arc of the teeth with a distance between 5-15 cm. See the positioning examples (p.X).
2. Plug the power adapter into the power network. Only one point/dot on the display lights up. Press START / STOP button and then the "Hi" (High) or "Lo" (Low) intensity indication appears.
3. The display shows "Hi" - a high intensity level. If you want to switch to the "Lo" (Low) intensity mode, press the TIMER (+) or TIMER (-) button once. Then the display will show "Lo". Every time before setting the required time, you will be able to select the Hi or Lo intensity level.
4. Press the START / STOP button once again and set time option will be displayed. Note that the device will save the time frame used in the latest bleaching procedure. Required time frame for operation is set by pressing TIMER + (increasing time) and TIMER- (reducing time):  
The display on the bleaching head will indicate timeframes between 10 sec. and 9 min. and 50 sec. with a comma between the two digits.  
Examples:  
Display indicates: 13. This means the timeframe is set to 13 min.  
Display indicates: 0.1. This means the timeframe is set to 10 sec.  
Display indicates: 9.5. This means the timeframe is set to 9 min. and 50 sec.  
The display on the bleaching head will indicate timeframes between 10 min. and 30 min. without a comma between the two digits.
5. Press the START / STOP button one final time and the device will start to emit light. A flashing dot on the display will indicate that the device is working properly.
6. If either of the TIMER + / TIMER - buttons is pressed during operation, the unit switches to Pause mode. A blinking display with the remaining time will indicate that the device is paused. Operation is resumed by pressing the TIMER + / TIMER button again.


7. To Stop the bleaching device at any time press the START / STOP button. If the device is not manually stopped it will stop after the selected time has elapsed. Note that the fans will continue to work for some time in order to cool down the unit.
8. At the end of the day turn off the device by holding down the Start / Stop button for 3-4 seconds or by disconnecting it from the power network.
9. In case of overheating the thermal protection is triggered and the device will cease to emit light.  
The display will indicate "Oh" and the fans will turn on for 1 minute.  
After cooling, the device will continue its normal operation.  
Overheating can only happen in an emergency situation or by accidental clogging of the cooling vents.  
That's why the cooling vents located on the bottom of the bleaching head must always be unobstructed in order to allow air circulation.


## VII. DAILY CARE AND MAINTENANCE

1. Disinfect the PVC window on the bleaching head before and after each patient with cotton or other nonabrasive cloth soaked in alcohol.
2. To disinfect the device and its parts, spray disinfectant onto a piece of soft cloth or cotton and clean it.  
Do not use abrasives or solvents as this may cause damage.

## POSSIBLE LAMP POSITIONING WITHIN DENTAL UNIT


**BLEACHING**


**STAND-BY**

**STAND-BY**

## ARCUS UNIT OPERATIVE POSITIONING



## VIII. PREPARATION FOR WORK WITH PATIENTS

1. Before starting work with the bleaching device, insulate the patient's non-calcified soft tissue, place a cover or protective napkin on the face, apply protective UV cream and put on protective glasses. The bleaching gel is applied according to the directions for use outlined by the gel manufacturer. The patient should be under constant supervision by the dentist during the bleaching process in order to avoid discomfort and to ensure all safety and technological procedures are applied.
2. It is desirable to check the sensitivity of the patient to the intensity of the light. The device has two levels of intensity (High and Low) and the optimal conditions are selected by adjusting the distance of the device to the patient's teeth. Testing should begin at a high level "Hi" (High) and a distance of 5 cm. If the thermal sensation is unpleasant and irritating, increase the distance to achieve an acceptable thermal comfort. If the thermal sensation is still unpleasant at 10-15 cm switch to low intensity- "Lo" (Low) and shorten the distance to the teeth. If the patient has no sensitivity or sensitivity cannot be measured for some reason, longer distance should be used.


## IX. POSITIONS OF BLEACHING HEAD AGAINST PATIENT'S TEETH

**CORRECT**


**INCORRECT**

## X. DETERMINATION OF DISTANCE TO PATIENT'S TEETH

### HIGH intensity

Hot

Warm


Acceptable

Weak

### LOW intensity

## XI. PROBLEMS AND SOLUTIONS

### The device does not work.

Check if the power cord from the adapter is plugged into the power network and all cable connectors are correctly coupled.

### Damaged power cord.

Do not use the device! Switch it off and immediately take it to the service facility

### Unpleasant or irritating thermal sensation by the patient.

Increase the distance of the device to the patient's teeth to 10-15 cm. If the thermal sensation is still unpleasant at this distance switch to low intensity- "Lo" (Low) and shorten the distance to the teeth.

### The display is blinking and no light is emitted.

The device is in PAUSE mode. To resume operation, press either + or - buttons on Bleaching Head.

### The Bleaching Head cannot be fixed.

Read Point 5. In Section V. INSTALLATION and slightly adjust the tension of the Bleaching Head until the desirable mobility is reached.

### The display only shows "Oh".

In case of overheating the thermal protection is triggered and the device will cease to emit light. The display will indicate "Oh" and the fans will turn on for 1 minute. After cooling, the device will continue its normal operation. Overheating can only happen in an emergency situation or by accidental clogging of the cooling vents. That's why the cooling vents located on the bottom of the bleaching head must always be unobstructed in order to allow air circulation.

For all other questions regarding the installation and/or operation of the device, please contact the manufacturer or your local authorized representative.

## XII. WARRANTY

1. The warranty period of the BLANCONE ARCUS BUILT-IN bleaching device is 24 (twenty-four) months from the date of purchase. If the purchase date is not filled in Section X. LED BLEACHUNG UNIT DATA, the warranty begins on the date of production.
2. During the warranty period, the replacement of the defective elements shall be done free of charge by the manufacturer. Note - LED modules have a 6-month warranty.
3. The device and its components must only be used as intended by the manufacturer purpose outlined in this User Manual. Any other use of the device will lead to a warranty lose and no obligations from the manufacturer.
4. If the device fails as a result of improper operation during the warranty period (mechanical, chemical, thermal, electrical), non-intended use, inappropriate storage etc..., warranty will be void and the repair cost will be at the expense of the user.

The device should not be used with damaged cables. If any such damages are noticed switch off the device and immediately take it to the service facility.

If liquids such as water and solvents, aggressive or flammable substances, vapors, insect or rodent enter the device, switch it off and immediately take it to the service facility. Damages from the above will void the warranty.

No claims shall be accepted for damages arising from electrical shocks, thunderstorms, non-compliance with electrotechnical safety measures or insufficient protection of patients, staff, other people, animals, plants, and objects from light radiation.

The manufacturer does not owe any compensation for lost profits during the period when the device is damaged or malfunctioned, no matter what the reason.

Claims for damages and claims arising out of non - compliance with the bleaching procedure are not accepted. Including but not limited to: shorter or longer time period use of the device with the bleaching material set by the gel manufacturer; unsatisfactory bleaching result; patient's harms from the bleaching material; overdosage; insufficient insulation of patient's non-calcified soft tissue; improper protection of the patient and staff; inappropriate, expired, intended for another wavelength and/or with unsuitable concentration bleaching material.

The warranty is forfeited and no claims shall be accepted for damages as a result of incorrect or insufficient care and attention for protection during transport, unpacking, moving, handling and storage of the device.

In the event of disputes arising out of the application and interpretation of this User Manual, those will be settled by the courts in the city of Plovdiv, under the current Bulgarian legislation.

5. The warranty of the device will be forfeited if any repairs or modifications are done by unauthorized personnel outside of the manufacturer's service facilities and/or non-original spare parts are used.
6. The manufacturer recommends that customers should check whether the device parameters are within the permissible limits once a year. Technical condition tests and verifications may only be carried out at the manufacturer's service facilities or by an authorized representative.
7. All repairs must be done at the manufacturer's service facilities with the following address:

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### XIII. SYMBOLS

#### SIMBOLI - SYMBOLE - SYMBOLES - SÍMBOLOS - SÍMBOLOS

	Manufacturer BG LIGHT LTD, 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria
	Pursuant to Directive 2012/19/EU, this symbol shows that the product must not be disposed of as urban waste at the end of its operating life.
	Pay attention to the instructions accompanied by this symbol.
	Fragile!
	Product with CE mark
	Applied part type B
$\Omega$	Ohm (electric resistance unit)
<b>s</b>	Second (time unit)
<b>W</b>	Watt (power unit)
<b>Hz</b>	Hertz (frequency unit)
<b>mm</b>	Millimeter (length unit)
<b>A</b>	Amper (electric current unit)
<b>V</b>	Volt (electric voltage unit)
<b>Pa</b>	Pascal (atmosphere pressure unit)
<b>°C</b>	Degrees Centigrade (temperature unit)
<b>g</b>	Gram (weight unit)
<b>SN</b>	Serial number
<b>LOT</b>	Lot

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Product: **BlancOne ARCUS\* Fixed**  
Classification: **Active device in Class IIa of the Medical Device Directive 93/42/EEC+AC:2007**

Classification is done by the manufacturer according to Medical Device Directive 93/42/EEC+AC:2007, Rule 9, Annex IX.

Notified body: TUV NORD Polska Sp. z o.o., ul.Mickiewicza 29, 40-085 Katowice, Poland.

This Declaration of conformity is valid only in combination with our certificates of Notified body TUV NORD Polska Sp. z o.o. Certificates N: AC090 100/0951/4047/2014, AC090 MD/0951/4047/2014, TNP/MDD/0122/4047/2015..

The manufacturer declares under its sole responsibility that the products are developed and produced in conformity with MDD 93/42/EEC+AC:2007 and the following applicable standards:

- EN 60601-1 Electrical safety,
- EN 60601-1-2 Electromagnetic compatibility
- EN 62304 Medical device software. Software life cycle processes.
- EN 62353 Medical electrical equipment. Recurrent test and test after repair of medical electrical equipment
- EN 62366 Medical devices. Application of usability engineering to medical devices.
- EN ISO 14155-1 Clinical investigation of medical devices for human subjects. General requirements.
- EN ISO 14155-2 Clinical investigation of medical devices for human subjects. Clinical investigation plans.
- EN ISO 14971 Medical devices – Application of risk management to medical devices
- EN ISO 15223-1 Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- EN 1041 Information supplied by the manufacturer with medical devices
- EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

All company products are manufactured under the current Quality Management System, ISO 9001 and ISO 13485.

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Date of production:

Date of purchase:

**The quick and gentle  
professional teeth whitening**

**[www.blancone.eu/pro](http://www.blancone.eu/pro)**

**IDS Spa**

17100 Savona - Via Valletta San Cristoforo, 28/10

Tel. +39 019 862080 - Fax +39 019 2304865

[www.ids dental.it](http://www.ids dental.it) - [info@ids dental.it](mailto:info@ids dental.it)

**BLANCONE®**

**The quick and gentle  
professional teeth whitening**

**[www.blancone.eu/pro](http://www.blancone.eu/pro)**

 **blancone**

 **BlancOne Team PRO**

 **blancone\_ids**

 **blancone**

**IDS Spa**

17100 Savona - Via Valletta San Cristoforo, 28/10  
Tel. +39 019 862080 - Fax +39 019 2304865  
[www.idsdental.it](http://www.idsdental.it) - [info@idsdental.it](mailto:info@idsdental.it)