



TePe Universal Care™

Handle	Polypropylene (PP).
Filaments	Polyamide (PA).
Neck	Can be angled after heating.
Autoclave	Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min, sterilization pressure: 3 bar / 300 kPa.
Pack	1 pc/blister pack, 12 blister packs/carton. The plastic front is made of Polyethylene terphthalate (PET). The backing card is made of cardboard. The carton is made of chlorine-free corrugated cardboard.
Recycling	The brush in the usual household waste. The blister pack should be separated and recycled: the plastic front in container for plastic waste and the backing card in container for recycling of cardboard. The carton should be placed in container for corrugated cardboard.
Product features	<ul style="list-style-type: none"> - This innovative brush with a unique angle provides easy access to the lingual and palatal surfaces of implant abutments. - Facilitates cleaning in very narrow areas. - Developed in collaboration with dental experts. - The angle of the neck can be adjusted; heat in hot water and bend.