



# TePe Denture Care

<b>Handle</b>	Polypropylene (PP).
<b>Filaments</b>	Polyamide (PA).
<b>Neck</b>	Can be angled after heating. Bend the handle backwards for improved grip.
<b>Pack</b>	1 pc/blister, 12 pcs/carton. Mixed colours. The carton is made of chlorine-free corrugated cardboard.
<b>Recycling</b>	The brush in the usual household waste. The blister pack should be separated and recycled: the plastic front in container for plastic waste and the backing card in container for recycling of cardboard. The carton should be placed in container for corrugated cardboard.
<b>Product features</b>	<ul style="list-style-type: none"><li>- Wavy trim for efficient cleaning.</li><li>- Designed as a conventional toothbrush.</li><li>- For improved grip, bend the handle backwards.</li></ul>