



# TePe Nova

<b>Handle</b>	Polypropylene (PP). Inlay: Thermoplastic elastomer (TPE).
<b>Filaments</b>	Polyamide (PA). The coloured filaments of the tip indicate the texture: Pink tip - x-soft, Yellow tip - soft, Blue tip - medium TePe Nova Medium has 1450 filaments/TePe Nova Soft has 1900 filaments/ TePe Nova Extra Soft has 2660 filaments.
<b>Neck</b>	Can be angled without heating.
<b>Autoclave</b>	Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min, sterilization pressure: 3 bar / 300 kPa.
<b>Pack</b>	1 pc/bag, 80 bags/carton. Mixed colours. The bag is made of Polypropylene (PP). 1 pc/blister pack, 12 blister packs/carton. Mixed colours. The plastic front is made of Polyethylene terphthalate (PET). The backing card is made of cardboard. The carton is made of chlorine-free corrugated cardboard.
<b>Recycling</b>	The toothbrush in the usual household waste and bag in container for plastic. The blister pack should be separated and recycled: the plastic front in container for plastic waste and the backing card in container for recycling of cardboard. The carton should be placed in container for corrugated cardboard.
<b>Product features</b>	<ul style="list-style-type: none"><li>- Tip for increased accessibility.</li><li>- End-rounded filaments for gentle cleaning.</li><li>- Tapered brush-head for improved access.</li><li>- Ergonomically designed non-slip handle for secure grip.</li></ul>