



TePe Special Care™

Handle	Polypropylene (PP).
Filaments	Extremely soft filaments in Polyamide (PA). TePe Special Care has 12000 filaments/TePe Special Care Compact has 7000 filaments.
Neck	Can be angled after heating.
Autoclave	Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min, sterilization pressure: 3 bar / 300 kPa.
Pack	1 pc/bag. 25 bags/carton. The bag is made of Polypropylene (PP). 1 pc/blister pack, 12 blister packs/carton. The plastic front is made of Polyethylene terphthalate (PET). The backing card is made of cardboard. The carton is made of chlorine-free corrugated cardboard.
Recycling	The brush and bag in the usual household waste. The blister pack should be separated and recycled: the plastic front in container for plasticwaste and the backing card in container for recycling of cardboard. The carton should be placed in container for corrugated cardboard.
Examples of areas of use	<ul style="list-style-type: none">- After oral surgery.- Extremely dry mouth.- Gingival lichen.- Canker sore and herpes.- Trauma.- Radiotherapy treatment.